

SC1005 Digital Logic

Simple Guide on Digital Oscilloscope

2023

Two models of oscilloscopes used in the lab

Tektronix TDS2012C



[user manual](#)

Rigol DS1104Z





[user manual](#)

We will use the oscilloscope to

- Display 2 waveforms
- Measure the voltage of a waveform
- Measure the frequency of a waveform

1. Turn on the main power and locate the On/Off button

Tektronix TDS2012C




Rigol DS1104Z




2. Press button to switch on the oscilloscope

Tektronix TDS2012C




Rigol DS1104Z



3. Connect signals to the oscilloscope

Tektronix TDS2012C



Channel 1 input

Channel 2 input

Rigol DS1104Z



Channel 1 input


Channel 2 input

4. Press AUTO button if no waveform is displayed on the oscilloscope

Tektronix TDS2012C




Rigol DS1104Z



5. Set probe to 1X for CH1 & CH2

Tektronix TDS2012C

(A) Press button to select Channel



(B) Press button to select Probe 1X

Rigol DS1104Z

(A) Press button to select Channel




(B) Press button to select Probe 1X

6. Read voltage levels – use squares on grid


Tektronix TDS2012C

Identify the ground level of each signal



Rigol DS1104Z


Identify the ground level of each signal



6. Read voltage levels – use squares on grid

Tektronix TDS2012C

Voltage = (no. of squares) x (volt per square)




CH1 voltage (approx.) = $1.6 \times 2 = 3.2\text{V}$

CH2 voltage (approx.) = $2.2 \times 2 = 4.4\text{V}$

Rigol DS1104Z

Voltage = (no. of squares) x (volt per square)




CH1 voltage (approx.) = $1.6 \times 2 = 3.2\text{V}$

CH2 voltage (approx.) = $2.2 \times 2 = 4.4\text{V}$

7. Read frequency – use squares on grid

Tektronix TDS2012C


Freq = $1/[(\text{no. of squares}) \times (\text{time per square})]$



CH1 & CH2 freq (approx.) = $1/(4 \times 250\text{us})$
= $1/1000\text{us} = 1 \text{ kHz}$

Rigol DS1104Z

Freq = $1/[(\text{no. of squares}) \times (\text{time per square})]$



CH1 & CH2 freq (approx.) = $1/(5 \times 200\text{us})$
= $1/1000\text{us} = 1 \text{ kHz}$

8. Read frequency – from oscilloscope

Tektronix TDS2012C

Freq = $1/[(\text{no. of squares}) \times (\text{time per square})]$



CH1 & CH2 freq (approx.) = $1/(4 \times 250\text{us})$
= $1/1000\text{us} = 1 \text{ kHz}$

Rigol DS1104Z


Freq = $1/[(\text{no. of squares}) \times (\text{time per square})]$



CH1 & CH2 freq (approx.) = $1/(5 \times 200\text{us})$
= $1/1000\text{us} = 1 \text{ kHz}$

9. Optional – adjust voltage scale


Tektronix TDS2012C



Adjust for Channel 1

Adjust for Channel 2

Rigol DS1104Z




Select Channel


Adjust for selected channel

10. Optional – adjust time scale

Tektronix TDS2012C




Rigol DS1104Z



11. Optional – adjust vertical positions of waveforms


Tektronix TDS2012C



Adjust for Channel 1

Adjust for Channel 2

Rigol DS1104Z



Select Channel

Adjust for selected channel

End of guide